

o r t h o p a e d i c s

VANCOGENX[®] Line



*Innovation,
Synergy,
Efficacy,
Evolution...
serving surgeons.*

TECRES[®] 
Advancing High Technology

www.tecres.it

Innovation


Tecres research laboratories developed the first **Vancomycin** and **Gentamicin** loaded bone cement.

VANCOGENX is the ideal solution for the fixation of antibiotic-loaded spacers and prosthetic components during revision surgery due to infections caused by resistant micro-organisms.

MRSA (Methicillin resistant *Staphylococcus aureus*) and MRSE (Methicillin resistant *Staphylococcus epidermidis*) infections, in particular, continue to present a serious, yet common problem.

VANCOGENX provides surgeons with the ideal solution for the treatment of these complex cases.

Proportion of methicillin resistant *staphylococcus auerus* (MRSA) isolates in participating countries in 2014.⁽¹⁴⁾



VANCOGENX features:

- allows effective, high-dose release of Vancomycin and Gentamicin
- guarantees excellent mechanical performances
- presents ideal viscosity for all types of use, for both manual and syringe application
- has a non-abrasive action due to the use of barium sulphate as the radiopaque tracer
- reduces the heat developed during curing and lower toxicity due to the low monomer content

Synergy


Vancomycin and Gentamicin

- Have a **synergic** action as their antibacterial efficacy is more potent when they are used together.^(1, 2, 3)
- Boast a spectrum of action covering approximately **90% of the pathogens** commonly isolated in orthopaedic infections.^(3, 4)
- Is **the most commonly used** and widely reported combination when bone cement is utilised in cases of concomitant treatment of infections.^(5, 6, 7)

Spectrum of action: microbiological efficacy

	Gram+					Gram-	
	Staphilococci MET-S	Staphilococci MET-R	Enterococci	Streptococci	Propionibacteria	Enterobacteria	<i>P. aeruginosa</i>
Gentamicin	good	low	good		no	high	
Vancomycin	high				good	no	
Vanco-Genta	synergic			high	good	synergic	

Pathogens isolated in prosthetic infections⁽⁴⁾



The advantages of VANCOGENX:


- broad spectrum of action
- synergic antimicrobial effect
- prevention of bacterial adhesion
- effective as a coadjuvant in the treatment of infection

Efficacy

VANCOGENX is effective:

in protecting the device

In vitro studies have shown that VANCOGENX exerts an anti-adhesion action by inhibiting bacterial roliferation.




The chart shows the adhesion of clinical isolates to PMMA loaded with Gentamicin, Vancomycin and a Vancomycin-Gentamicin combination (VANCOGENX).

Use of the Vanco-Genta combination prevents bacteria from adhering to the PMMA⁽³⁾


as a coadjuvant in the treatment of infection

*In vivo** studies have proven VANCOGENX's therapeutic efficacy in experimental osteomyelitis models.

Microbiological analysis of culture smears



Bacterial load in the bone



Description of the study:

following induction of MRSA infection for 4 weeks, animals underwent debridement and were assigned to 4 groups representing different treatment options: stainless steel nail (group 1); Vancogenx-coated stainless steel nail (group 2); no treatment (group 3); systemic treatment with I.V. Teicoplanin. After 3 weeks, smears and bone fragment specimens were harvested to evaluate the progression of the infection and the effect of/on bacterial load was evaluated following exposure to the various treatment options using both culture swabs and bone analysis.

Local therapy with Vancogenx eradicates infection and provides better results than systemic therapy with teicoplanin ⁽¹⁵⁾.

* *in vivo* animal studies

Evolution

VANCOGENX-SPACE devices are the exclusive preformed spacers that bring together the extraordinary mechanical and pharmacological characteristics of Tecres spacers and the efficacy of the Vancomycin-Gentamicin combination. They are indicated for temporary artificial joint replacement when prosthetic joints require removal due to sepsis.

The microbiological efficacy of VANCOGENX-SPACE devices can be attributed to the special production process developed by Tecres that increases the device's porosity without compromising its mechanical performance.⁽⁸⁾

This permits antibiotic release that is:


- sustained and high-dose
- calibrated
- superior to that achieved with theatre-made devices using antibiotic-loaded bone cement.



VANCOGENX-SPACE features:

- **boast a spectrum of action** that covers approximately 90% of the pathogens commonly isolated in orthopaedic infections.
- **prevent bacterial colonisation** and proliferation by guaranteeing the continuous presence of antibiotic agents on surface of the device ⁽³⁾
- **reduce the risk of bacterial resistance.** They permit a consistent initial release of antibiotic that is kept at a high level and sustained over time


In vitro microbiological release (V+G)



Concentrations above the minimum inhibitory concentration (MIC) for the most common pathogens are maintained for at least 30 days (MIC > 10µg/ml). ⁽¹⁰⁾


VANCOGENX-SPACE Hip

Comes in 6 sizes. (3 head sizes with standard stem and 3 head sizes with XL stem).
If the stem is unstable and when using Vancogenx Space Hip XL, proximal cementing of the neck should be performed using Vancogenx cement.



VANCOGENX-SPACE KNEE

Comes in 4 sizes. Requires component cementing using Vancogenx cement.



Ordering information

Spacer for HIP

Spacer G - Gentamicin Loaded			
Stem	Head Size (mm)	Round Stem	Flat Stem
Short Stem	46	SPC46/G	SPC0620
	54	SPC54/G	SPC0720
	60	SPC60/G	SPC0820
Long Stem	46	SPC46/GXL	SPC0920
	54	SPC54/GXL	SPC1020
	60	SPC60/GXL	SPC1120

Vancogenx Space Hip - Genta + Vanco Loaded			
Stem	Head Size (mm)	Round Stem	Flat Stem
Short Stem	46	SPC0030	SPC0630
	54	SPC0130	SPC0730
	60	SPC0230	SPC0830
Long Stem	46	SPC0330	SPC0930
	54	SPC0430	SPC1030
	60	SPC0530	SPC1130

Trial Set (3-size set)	
Short Round Stem	SPG03
Long Round Stem	SPG03XL
Short Flat Stem	SPC90Z0
Long Flat Stem	SPC91Z0

VANCOGENX

Code	Product	Details
12A2520	Vancogenx	40 g

Spacer for KNEE

Spacer K - Gentamicin Loaded	
Tibial Dimension (mm)	Code
60 – small	SPK6054/G
70 – medium	SPK7064/G
80 – large	SPK8074/G

Vancogenx Space K - Genta + Vanco Loaded	
Tibial Dimension (mm)	Code
60 – small	SPK0030
70 – medium	SPK0130
80 – large	SPK0230
90 – extralarge	SPK0330

Trial Set	
Trial Set (small–medium large)	SPK03
Trial extralarge	SPK03Z0

Spacer for SHOULDER

Spacer S - Gentamicin Loaded	
Head Size (mm)	Code
41	SPS0020
46	SPS46/G

Trial	
Two–size	SPS90Z0

Bibliography

1. Watanakunakorn C, Bakie C.
Synergism of vancomycin-gentamicin and vancomycin-streptomycin against enterococci.
Antimicrob Agents Chemother. 1973 Aug;4(2):120-4.
2. Watanakunakorn C, Tisone JC.
Synergism between vancomycin and gentamicin or tobramycin for methicillin-susceptible and methicillin-resistant *Staphylococcus aureus* strains.
Antimicrob Agents Chemother. 1982 Nov;22(5):903-5.
3. E. Bertazzoni Minelli, T. Della Bora, A. Benini
Different microbial biofilm formation on polymethylmethacrylate (PMMA) bone cement loaded with gentamicin and vancomycin
Anaerobe. 2011 Dec; 17(6): 380-3
4. Trampuz A, Zimmerli W.
Prosthetic joint infections: update in diagnosis and treatment.
Swiss Med Wkly. 2005 Apr 30;135(17-18):243-51. Review.
5. Penner MJ, Masri BA, Duncan CP.
Elution characteristics of vancomycin and tobramycin combined in acrylic bone-cement.
J Arthroplasty. 1996 Dec;11(8):939-44.
6. Masri BA, Duncan CP, Beauchamp CP.
Long-term elution of antibiotics from bone-cement: an in vivo study using the prosthesis of antibiotic-loaded acrylic cement (PROSTALAC) system.
J Arthroplasty. 1998 Apr;13(3):331-8.
7. Bertazzoni Minelli E, Caveiari C, Benini A.
Release of antibiotics from polymethylmethacrylate cement.
J Chemother. 2002 Oct;14(5):492-500.
8. Soffiatti R.
In "Infection and local treatment in orthopedic surgery"
2007, Springer Verlag. *The preformed spacers from the idea to the realization of industrial device.*
9. Springer BD et al.
Systemic safety of high-dose antibiotic-loaded cement spacers after resection of an infected total knee arthroplasty.
Clin Orthop Relat Res. 2004 Oct;(427):47-51.
10. Bertazzoni Minelli E et al.
PK-PD parameters ratio of antibiotics released from antibiotic-loaded spacers in drainage fluids
ESCMID Munich, Germany 2007.
11. Bertazzoni Minelli E., Benini A.
In "Infection and local treatment in orthopedic surgery"
2007, Springer Verlag. *The Gentamicin.-Vancomycin Spacer: A Pharmacological study.*
12. Garvin KL, Hanssen AD.
Infection after total hip arthroplasty. Past, present, and future.
J Bone Joint Surg Am. 1995 Oct;77(10):1576-88. Review.
13. Moojen DJ, Hentenaar B, Charles Vogely H, Verbout AJ, Castelein RM, Dhert WJ.
In vitro release of antibiotics from commercial PMMA beads and articulating hip spacers.
J Arthroplasty. 2008 Dec;23(8):1152-6. Epub 2008 Mar 4.
14. EARSS Annual Report 2014.
Antimicrobial resistance surveillance in Europe.
15. Giavaresi G. et al.
Preliminary investigations on a new Genta and Vanco-coated PMMA nail for the treatment of bone and intramedullary infections: an experimental study in the rabbit.
J Orthop Res. 2008. n; 26(6):785-92.

