


Tips when using ECG247



- Download the ECG 247 app and follow the instructions carefully
- Do not remove the green foil over the battery of the electrode patch
- Hair must be removed and the skin over the sternum must be clean and dry
- The sensor is glued in the middle of the sternum a little below the throat pits
- Keep your phone nearby for the entire test period
- The sensor usually automatically reconnects to the phone if it has been out of range
- Turn off battery saving features and charge your phone every day

Further advice:

- You should treat the sensor carefully during the test period so that it does not come loose
- Check that the test runs as normal in the ECG247 app at least twice a day. If the test has stopped, you can start a new test in the app without removing the sensor
- You can exercise and shower with the sensor on, but you cannot bathe/swim
- By pressing briefly and hard in the middle of the sensor, a manually ECG recording is stored
- You stop the test in the ECG247 app. Recommended test duration is at least 3 days
- The test can be shared with your doctor by pressing the share button in the app.