



Create the health with science and technology

Application

- Diagnose the pregnancy of Swine, Bovine, Ovine, Canine, Feline and other animal
- Diagnose animal's uteri and fetus condition
- Automatically calculate the Back fat and Lean percentage

Function feature

- Display mode: B、2B、4B、B/M、M, total 5 kinds
- 100 frames permanent storage, 256 frames cine loop
- Internal eight pseudo colors (incl.B/W)
- OB measurement: EDD and GA for Bovine, equine, ovine, canine, feline, goat, swine and llama
- OB report
- Grid: 6 kinds (None Grid、Line grid、Point grid、10mm spcae point ruler、2mm spcae point ruler、10mm and 2mm spcae point ruler)
- USB port support large capacity USB disk (support document management /software upgrades/one-key storage)
- Support color laser printer directly print various image and reports


Standard configuration

- Main unit: 1 pc
- Probe : 1 pc
- Lithium battery : 1 pc
- Sunshine cover : 1 unit
- Black plastic suitcase : 1 pc

Optional configuration

- Lithium battery
- Leather Case saddlebag

Optional probe



6.5MHz multi-frequency linear rectal probe



4.0MHz multi-frequency convex rectal probe (patent probe)



3.5MHz Veterinary loin muscle linear probe



3.5MHz multi-frequency convex probe



7.5MHz multi-frequency linear probe



6.5MHz multi-frequency micro-convex probe



RKU10

Veterinary B mode ultrasound scanner



KAI XIN

XUZHOU KAIXIN ELECTRONIC INSTRUMENT CO.,LTD.

Add: Kaixin Mansion, C-01, Economic Development Zone, Xuzhou, Jiangsu, China.

Tel: +86-516-87732932/87733758 Fax: +86-516-87732932/87792848 Zip Code: 221004

Http://www.kxele.com Business: info@kxele.com Service: service@kxele.com

DISTRIBUTOR

THE PRODUCT TAKES REAL OBJECT AS QUASI, IF ANY FLUCTUATION HAPPENS, FORGIVING NO ANOTHER NOTICE IN THE LINE.

Edi tion No: ZL-RKU10-2013-01

KAI XIN

Professional manufacturer of ultrasound scanner for veterinary use



RKU10

Veterinary B mode ultrasound scanner



With Most Convenient Experience You Never Know!

RKU adopt sophisticated hardware craft, ultra light weight, ultra long-time battery work, unique ring design, multiple use method, convenient for long-time using but eliminating operate fatigue.



Sophisticated craft, have soft rubber coat, splash-proof, dust-proof. Easy to clean



Multiple use method: hand-held, chest hanging, ect.



Available for left/right hand alternative use, eliminating operate fatigue.



Net weight: 1.1kg. Smart and light.



5.7" High resolution color LED, high brightness, high contrast, wide visual.

Features



High-capacity lithium Battery, more than 6 hours work.