



○ : FDA Approved X : Not Approved △ : Not Tested

	Disinfectant							
	Cidex OPA	Gigasept FF	Virkon	Cidex Plus	Wavicide-01	AIDAL PLUS	Cetylclde-G	Sporicidin
C1-6	○	○	○	○	○	○	○	○
C1-6T	○	X	○	○	△	○	△	X
C1-6i	○	○	○	○	○	○	○	○
SC1-6	○	○	○	○	○	○	○	○
C5-8	○	○	△	○	○	△	△	△
C5-8N	○	X	○	○	△	○	△	X
SC1-4HS	○	○	○	○	○	○	○	○
SC1-4H	○	○	○	○	○	○	○	○
SC1-6H	○	○	○	○	○	○	○	○
L3-8	○	○	○	○	○	○	○	○
L3-8i	○	○	○	○	○	○	○	○
L3-12	○	○	○	○	○	○	○	○
L3-12T	○	X	○	○	△	○	△	X
L3-12i	○	○	○	○	○	○	○	○
L3-12H	○	○	○	○	○	○	○	○
L3-12H ^{WD}	○	○	○	○	○	○	○	○
L3-12X	○	○	○	○	○	○	○	○
L8-17	○	○	○	○	○	○	○	○
L8-17i	○	○	○	○	○	○	○	○
L8-17X	○	○	○	○	○	○	○	○
MP1-5X	○	○	○	○	○	○	○	○
SP1-5	○	○	○	○	○	○	○	○
SP1-5i	○	○	○	○	○	○	○	○
SP3-8	○	○	○	○	○	○	○	○
SP1-5X	○	○	○	○	○	○	○	○
VC1-6	○	○	○	○	△	△	△	△
SVC1-6	○	○	○	○	△	△	△	△
SVC1-6H	○	○	○	○	△	△	△	△
E3-10	○	△	○	○	△	△	△	△
EN3-10	○	△	△	○	○	△	△	△
EV3-10	○	○	○	○	○	○	○	○
EV3-10T	○	○	○	○	○	○	○	○
EC3-10	○	○	○	○	○	○	○	○
EC3-10T	○	○	○	○	○	○	○	○
E3-10H	○	△	○	○	△	△	△	△
EC3-10H	○	○	○	○	○	○	○	○
EV3-10H	○	○	○	○	○	○	○	○
IO3-12	○	○	○	○	○	○	○	○
CW2.0	△	△	△	△	△	△	△	△
CW5.0	△	△	△	△	△	△	△	△



○ : FDA Approved X : Not Approved △ : Not Tested

	Pre-Cleaner (1/2)									
	T-spray	T-spray II	Enzol	Theracide Plus(S5N)	Transeptic	Septihol	Metrizyme	Klenzyme	McKesson	Cidezime
C1-6	△	○	○	○	○	○	△	△	△	△
C1-6T	△	△	○	△	△	△	△	△	△	△
C1-6i	△	○	○	○	○	○	△	△	△	△
SC1-6	△	○	○	○	○	○	△	△	△	△
C5-8	○	○	△	△	○	△	△	△	△	△
C5-8N	△	△	○	△	△	△	△	△	△	△
SC1-4HS	△	○	○	○	○	○	△	△	△	△
SC1-4H	△	○	○	○	○	○	△	△	△	△
SC1-6H	△	○	○	○	○	○	△	△	△	△
L3-8	△	○	○	○	○	○	△	△	△	△
L3-8i	△	○	○	○	○	○	△	△	△	△
L3-12	△	○	○	○	○	○	△	△	△	△
L3-12T	△	△	○	△	△	△	△	△	△	△
L3-12i	△	○	○	○	○	○	△	△	△	△
L3-12H	△	○	○	○	○	○	△	△	△	△
L3-12H ^{WD}	△	○	○	○	○	○	△	△	△	△
L3-12X	△	○	○	○	○	○	△	△	△	△
L8-17	△	○	○	○	○	○	△	△	△	△
L8-17i	△	○	○	○	○	○	△	△	△	△
L8-17X	△	○	○	○	○	○	△	△	△	△
MP1-5X	△	○	○	○	○	○	△	△	△	△
SP1-5	△	○	○	○	○	○	△	△	△	△
SP1-5i	△	○	○	○	○	○	△	△	△	△
SP3-8	△	○	○	○	○	○	△	△	△	△
SP1-5X	△	○	○	○	○	○	△	△	△	△
VC1-6	△	○	○	○	○	△	△	△	△	△
SVC1-6	△	○	○	○	○	△	△	△	△	△
SVC1-6H	△	○	○	○	○	△	△	△	△	△
E3-10	○	△	△	△	○	△	△	△	△	△
EN3-10	○	△	△	△	△	△	○	○	○	○
EV3-10	△	○	○	○	○	○	△	△	△	△
EV3-10T	△	○	○	○	○	○	△	△	△	△
EC3-10	△	○	○	○	○	○	△	△	△	△
EC3-10T	△	○	○	○	○	○	△	△	△	△
E3-10H	○	△	△	△	○	△	△	△	△	△
EC3-10H	△	○	○	○	○	○	△	△	△	△
EV3-10H	△	○	○	○	○	○	△	△	△	△
IO3-12	△	○	○	○	○	○	△	△	△	△
CW2.0	△	△	△	△	△	△	△	△	△	△
CW5.0	△	△	△	△	△	△	△	△	△	△



○ : FDA Approved X : Not Approved △ : Not Tested

	Pre-Cleaner (2/2)			Gel						
	Endozone	Super Sani-Cloth	Cleanisep wipes	Aquasonic 100	Soflean Gel	Care Sonic Gel	Ultrasound gel	Ultrasound Transmission gel	Scan	Sonogel
C1-6	○	○	○	○	○	○	○	△	△	△
C1-6T	△	○	○	△	△	△	○	○	△	△
C1-6i	○	○	○	○	○	○	○	△	△	△
SC1-6	○	○	○	○	○	○	○	△	△	△
C5-8	△	△	△	○	△	△	△	△	○	○
C5-8N	△	○	○	△	△	△	○	○	△	△
SC1-4HS	○	○	○	○	○	○	○	△	△	△
SC1-4H	○	○	○	○	○	○	○	△	△	△
SC1-6H	○	○	○	○	○	○	○	△	△	△
L3-8	○	○	○	○	○	○	○	△	△	△
L3-8i	○	○	○	○	○	○	○	△	△	△
L3-12	○	○	○	○	○	○	○	△	△	△
L3-12T	△	○	○	△	△	△	○	○	△	△
L3-12i	○	○	○	○	○	○	○	△	△	△
L3-12H	○	○	○	○	○	○	○	△	△	△
L3-12H ^{WD}	○	○	○	○	○	○	○	△	△	△
L3-12X	○	○	○	○	○	○	○	△	△	△
L8-17	○	○	○	○	○	○	○	△	△	△
L8-17i	○	○	○	○	○	○	○	△	△	△
L8-17X	○	○	○	○	○	○	○	△	△	△
MP1-5X	○	○	○	○	○	○	○	△	△	△
SP1-5	○	○	○	○	○	○	○	△	△	△
SP1-5i	○	○	○	○	○	○	○	△	△	△
SP3-8	○	○	○	○	○	○	○	△	△	△
SP1-5X	○	○	○	○	○	○	○	△	△	△
VC1-6	△	○	△	△	△	△	△	△	△	△
SVC1-6	△	○	△	△	△	△	△	△	△	△
SVC1-6H	△	○	△	△	△	△	△	△	△	△
E3-10	△	○	△	○	△	△	△	△	△	△
EN3-10	△	○	△	△	△	△	△	△	△	△
EV3-10	○	○	○	○	○	○	○	△	△	△
EV3-10T	○	○	○	○	○	○	○	△	△	△
EC3-10	○	○	○	○	○	○	○	△	△	△
EC3-10T	○	○	○	○	○	○	○	△	△	△
E3-10H	△	○	△	○	△	△	△	△	△	△
EC3-10H	○	○	○	○	○	○	○	△	△	△
EV3-10H	○	○	○	○	○	○	○	△	△	△
IO3-12	○	○	○	○	○	○	○	△	△	△
CW2.0	△	△	△	△	△	△	△	△	△	△
CW5.0	△	△	△	△	△	△	△	△	△	△